

Anti-CD51 antibody [23C6] (STJ16100100)

STJ16100100

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD51 is suitable for use in Immunofluorescence, Immunoprecipitation, Flow Cytometry and Immunohistochemistry research applications.
Applications	IF/IP/FC/IHC-F
Host/Source	Mouse
Reactivity	Human/Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	23C6
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:40-1:80
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Osteoclasts from osteoclastomas.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	