

Anti-CD15 antibody [V1MC6] {FITC} (STJ16100109)
STJ16100109**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD15 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	V1MC6
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 4mg/ml BSA and 0.1% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human CD15
Immunogen Region	
Specificity	
Immunogen Sequence	