

Anti-CD64 antibody [10.1] {PE} (STJ16100120)
STJ16100120**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD64 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10.1
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with BSA, sucrose and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human CD64
Immunogen Region	
Specificity	
Immunogen Sequence	