

Anti-Cdk4 antibody [DCS-31] (STJ16100131)

STJ16100131

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cdk4 is suitable for use in Immunocytochemistry, Immunofluorescence, Western Blot, Flow Cytometry and Immunoprecipitation research applications.
Applications	ICC/IF/WB/FC/IP
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DCS-31
Concentration	0.2 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.2% BSA and 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	1019
Gene Symbol	CDK4
Uniprot ID	CDK4_HUMAN
Immunogen	Purified recombinant cdk4 protein
Immunogen Region	
Specificity	
Immunogen Sequence	