

Anti-Nucleophosmin antibody [FC-61991] (STJ16100133)

STJ16100133

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Nucleophosmin is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry, Immunofluorescence, ELISA and Immunoprecipitation research applications.
Applications	WB/IHC/ICC/IF/ELISA/IP
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FC-61991
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA 1-2 µg/mL ICC 2-10 µg/mLIF 2-10 µg/mL
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4869
Gene Symbol	NPM1
Uniprot ID	NPM_HUMAN
Immunogen	NPM/B23 purified from rat hepatoma
Immunogen Region	
Specificity	
Immunogen Sequence	