

Anti-CD46 antibody [MEM-258] (STJ16100134)

STJ16100134

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD46 is suitable for use in Flow Cytometry, Immunoprecipitation and Western Blot research applications.
Applications	FC/IP/WB
Host/Source	Mouse
Reactivity	Cow/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-258
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4179
Gene Symbol	CD46
Uniprot ID	MCP_HUMAN
Immunogen	HPB-ALL human T cell line
Immunogen Region	
Specificity	
Immunogen	
Sequence	