

Anti-CD46 antibody [MEM-258] (STJ16100135)

STJ16100135

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD46 is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P/IHC-F/WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-258
Concentration	5 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:10 IHC-F 1:10 WB ELISA
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4179
Gene Symbol	CD46
Uniprot ID	MCP_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

