

Anti-CD95 antibody [B-R18] (STJ16100141)

STJ16100141

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD95 is suitable for use in Flow Cytometry, Immunofluorescence and Immunohistochemistry research applications.
Applications	FC/IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B-R18
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution	Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30 min at
Range	RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	356
Gene Symbol	FASLG
Uniprot ID	TNFL6_HUMAN
Immunogen	A BALB/c mouse was immunized with recombinant Fas antigen. Fusion partner: P3-X63-Ag8.653.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	