

Anti-Lipid A antibody [26-5] (STJ16100142)

STJ16100142

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lipid A is suitable for use in ELISA and AGG research applications.
Applications	ELISA/AGG
Host/Source	Mouse
Reactivity	Bacteria

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	26-5
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide and 0.1% BSA
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic Lipid A
Immunogen Region	
Specificity	
Immunogen Sequence	