

Anti-CD100 antibody [133.1G6] (STJ16100144)

STJ16100144

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD100 is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC-F
Host/Source	Mouse
Reactivity	Human/Cow/Baboon

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	133.1G6
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	10507
Gene Symbol	SEMA4D
Uniprot ID	SEM4D_HUMAN
Immunogen	PHA stimulated PBL
Immunogen Region	
Specificity	
Immunogen Sequence	