

Anti-H-FABP antibody [6+60E+03] (STJ16100146)

STJ16100146

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-H-FABP is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat/Swine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6+60E+03
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	2170
Gene Symbol	FABP3
Uniprot ID	FABPH_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	