

Anti-Alpha-v Beta Integrin antibody [BV3] (STJ16100153)

STJ16100153

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha-v Beta Integrin is suitable for use in Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation and Immunohistochemistry research applications.
Applications	FC/ELISA/IF/IP/IHC-P
Host/Source	Mouse
Reactivity	Human/Chicken/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV3
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 ELISA 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence