

Anti-CD4 antibody [MEM-241] {FITC} (STJ16100167)

STJ16100167

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-241
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
Uniprot ID	CD4_HUMAN
Immunogen	2 N-terminal domains of human CD4 fused to human IgG1 Fc
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081