

Anti-Phospho-Tyrosine antibody [P-TYR-01] (STJ16100174)

STJ16100174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Phospho-Tyrosine is suitable for use in Western Blot, Immunocytochemistry and Flow Cytometry research applications.
Applications	WB/ICC/FC
Host/Source	Mouse
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	P-TYR-01
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Phosphotyrosine conjugated to bovine serum albumin.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	