

Anti-CD177 antibody [MEM-166] (STJ16100178)

STJ16100178

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD177 is suitable for use in Flow Cytometry, Immunoprecipitation and Western Blot research applications.
Applications	FC/IP/WB
Host/Source	Mouse
Reactivity	Non-Human Primates/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-166
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	57126
Gene Symbol	CD177
Uniprot ID	CD177_HUMAN
Immunogen	Human granulocytes
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081