

Anti-Proteinase Inhibitor 9 antibody [PI9-17] (STJ16100180)

STJ16100180

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Proteinase Inhibitor 9 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PI9-17
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:100 WB
Formulation	serumfree culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	fl recombinant PI9
Immunogen Region	
Specificity	
Immunogen Sequence	