

Anti-CD44 antibody [MEM-263] (STJ16100181)

STJ16100181

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD44 is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	FC/IP/WB/IHC-P
Host/Source	Mouse
Reactivity	Pig/Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-263
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	COS-7 cells (African Green Monkey).
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081