

Anti-CD44 antibody [MEM-263] {PE} (STJ16100182)

STJ16100182

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD44 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Pig/Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-263
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	COS-7 cells (African Green Monkey).
Immunogen Region	
Specificity	
Immunogen Sequence	