

Anti-CD80 antibody [MEM-233] (STJ16100188)

STJ16100188

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD80 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-233
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	941
Gene Symbol	CD80
Uniprot ID	CD80_HUMAN
Immunogen	Extracellular domain of human CD80 fused to human IgG1 (Fc)
Immunogen Region	
Specificity	
Immunogen	
Sequence	