

Anti-HLA-G antibody [MEM-G/9] {FITC} (STJ16100198)

STJ16100198

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-G is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-G/9
Concentration	1 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3135
Gene Symbol	HLA-G
Uniprot ID	HLAG_HUMAN
Immunogen	Recombinant human HLA-G refolded with beta2-microglobulin and peptide.
Immunogen Region	
Specificity	
Immunogen Sequence	