

Anti-HLA-E antibody [MEM-E/02] (STJ16100201)

STJ16100201

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-HLA-E is suitable for use in Western Blot and Immunohistochemistry research applications. |
| Applications | WB/IHC-P |
| Host/Source | Mouse |
| Reactivity | Human/Non-Human Primates |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | MEM-E/02 |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|-------------------------|--|
| Gene ID | 3133 |
| Gene Symbol | HLA-E |
| Uniprot ID | HLAE_HUMAN |
| Immunogen | Recombinant human HLA-E denaturated heavy chain. |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |