

Anti-HLA-Class I antibody [MEM-123] (STJ16100203)

STJ16100203

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-Class I is suitable for use in Flow Cytometry, Immunoprecipitation and ELISA research applications.
Applications	FC/IP/ELISA
Host/Source	Mouse
Reactivity	Cow/Human/Non-Human Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-123
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	COS-7 African green monkey kidney cells
Immunogen Region	
Specificity	
Immunogen	
Sequence	