

Anti-NTAL antibody [NAP-03] (STJ16100206)

STJ16100206

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-NTAL is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	IP/WB/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NAP-03
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7462
Gene Symbol	LAT2
Uniprot ID	NTAL_HUMAN
Immunogen	Peptide corresponding to aa 196-212 of human NTAL coupled to KLH.
Immunogen Region	
Specificity	
Immunogen Sequence	