

Anti-STAT 1 antibody [SM1] (STJ16100213)

STJ16100213

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-STAT 1 is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	IP/WB/IHC-P
Host/Source	Mouse
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM1
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	STAT1 peptide sequence 721-733 (DNLLPMSPEEFDE).
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

