

Anti-Nectin-2 antibody [502-57] (STJ16100216)

STJ16100216

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Nectin-2 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IP/WB
Host/Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	502-57
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IHC 1:50 ELISA 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	5819
Gene Symbol	NECTIN2
Uniprot ID	NECT2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	