

Anti-TLR9 antibody [5G5] (STJ16100217)

STJ16100217

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR9 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/IF/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5G5
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IHC 1:50 WB 1:50 IF 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	54106
Gene Symbol	TLR9
Uniprot ID	TLR9_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	