

## Anti-Smoothelin antibody [R4A] (STJ16100219)

STJ16100219

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Smoothelin is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	ICC/IHC-F/IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Dog/Chicken/Feline/Monkey/Swine

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	R4A
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	6525
<b>Gene Symbol</b>	SMTN
<b>Uniprot ID</b>	SMTN_HUMAN
<b>Immunogen</b>	cytoskeletal extract of chicken gizzard
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	

