

## Anti-Pan-Keratin antibody [80] (STJ16100232)

STJ16100232

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Pan-Keratin is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
<b>Applications</b>	IHC-F/WB/FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	80
<b>Concentration</b>	15 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC 1:50 IHC 1:10-1:50
<b>Formulation</b>	Culture medium with 0.01M Sodium Azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Callus cyokeratins from human skin
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	

