

Anti-Vimentin antibody [V9] (STJ16100234)

STJ16100234

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vimentin is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC-P/IHC-F/WB/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	V9
Concentration	25 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:300 IHC-P 1:50-1:100
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
Uniprot ID	VIME_HUMAN
Immunogen	Vimentin isolated from cow lens
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081