

Anti-Keratin 5+8 antibody [RCK102] (STJ16100236)

STJ16100236

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 5+8 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/IHC-P/IF/WB
Host/Source	Mouse
Reactivity	Dog/Hamster/Human/Mouse/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RCK102
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Lung cancer cell preparation
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	