

Anti-Neurofilament antibody [NF-01] (STJ16100237)

STJ16100237

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neurofilament is suitable for use in Immunohistochemistry, Western Blot and Immunocytochemistry research applications.
Applications	IHC-P/WB/ICC
Host/Source	Mouse
Reactivity	Mammalian

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NF-01
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081