

Anti-Keratin 13 antibody [2D7] (STJ16100243)

STJ16100243

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 13 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2D7
Concentration	15 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:5-1:25
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3860
Gene Symbol	KRT13
Uniprot ID	K1C13_HUMAN
Immunogen	Cytokeratin purified from human esophagus
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081