

Anti-Keratin 10 antibody [DE-K10] (STJ16100249)

STJ16100249

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 10 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DE-K10
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:100-1:250 IHC 1:5-1:20
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3858
Gene Symbol	KRT10
Uniprot ID	K1C10_HUMAN
Immunogen	Cytoskeletal preparation from human epidermis
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081