

Anti-Keratin 18 antibody [RGE 53] (STJ16100250)

STJ16100250

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 18 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/WB
Host/Source	Mouse
Reactivity	Chicken/Dog/Human/Mouse/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RGE 53
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:5-1:10
Formulation	Contains 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3875
Gene Symbol	KRT18
Uniprot ID	K1C18_HUMAN
Immunogen	Isolated from HeLa cells
Immunogen Region	
Specificity	
Immunogen	
Sequence	