

Anti-Keratin 18 antibody [RCK106] (STJ16100251)

STJ16100251

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 18 is suitable for use in Flow Cytometry, Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	FC/ICC/IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RCK106
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3875
Gene Symbol	KRT18
Uniprot ID	K1C18_HUMAN
Immunogen	Cytokeratins from the Human bladder carcinoma cell line T24
Immunogen Region	
Specificity	
Immunogen Sequence	

