

## Anti-Cardiotin antibody [SR-2] (STJ16100265)

STJ16100265

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cardiotin is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-F/IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Swine

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SR-2
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-F 1:25-1:100 IHC-P 1:25-1:100 WB 1:50-1:500
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgM
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	100 kDa cardiotin subunit
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

