

Anti-JAM-1 antibody [BV16] (STJ16100272)

STJ16100272

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-JAM-1 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/FC/IF/IP/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV16
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IP 1:50 IHC 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
Uniprot ID	JAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	