

**Anti-JAM-1 antibody [M.Ab.F11] (STJ16100273)**

STJ16100273

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-JAM-1 is suitable for use in Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	FC/ELISA/IF/IP/IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M.Ab.F11
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1:50 IP 1:50 IHC 1:50 WB 1:50 IP 1:50
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">50848</a>
<b>Gene Symbol</b>	<a href="#">F11R</a>
<b>Uniprot ID</b>	<a href="#">JAM1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	