

Anti-TLR3 antibody [TLR3.7] (STJ16100285)

STJ16100285

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR3 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunofluorescence, Immunoprecipitation and Immunohistochemistry research applications.
Applications	IHC-F/FC/IF/IP/IHC-P
Host/Source	Mouse
Reactivity	Human/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TLR3.7
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IP 1:50 IHC 1:50 IP 1:50
Formulation	PBS with 0.1 % BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	