

Anti-Keratin 6 antibody [Ks6.KA12] (STJ16100287)

STJ16100287

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 6 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Ks6.KA12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:10
Formulation	PBS with 0.5% BSA and 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Keratin K6 of human callus cytoskeletal preparation
Immunogen Region	
Specificity	
Immunogen Sequence	