

Anti-Macrophages antibody [25F9] (STJ16100303)

STJ16100303

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Macrophages is suitable for use in Immunohistochemistry, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	25F9
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:10 WB 1:10 IHC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence