

Anti-F4/80 antibody [BM8] (STJ16100305)

STJ16100305

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-F4/80 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/IHC-P/WB
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BM8
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:50 FC 1:50 IHC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence