

**Anti-CD61 antibody [2f2] (STJ16100340)**  
STJ16100340**GENERAL INFORMATION**


|                          |                                                                                                        |
|--------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                     |
| <b>Short Description</b> | Mouse monoclonal antibody anti-CD61 is suitable for use in Immunohistochemistry research applications. |
| <b>Applications</b>      | IHC-P                                                                                                  |
| <b>Host/Source</b>       | Mouse                                                                                                  |
| <b>Reactivity</b>        | Human                                                                                                  |

**PRODUCT PROPERTIES**

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>Clone ID</b>            | 2f2                                                               |
| <b>Concentration</b>       |                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Purification</b>        |                                                                   |
| <b>Dilution Range</b>      | IHC 1: 25 – 1:100                                                 |
| <b>Formulation</b>         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| <b>Isotype</b>             | IgG1k                                                             |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt.                                      |

**TARGET INFORMATION**

|                           |
|---------------------------|
| <b>Gene ID</b>            |
| <b>Gene Symbol</b>        |
| <b>Uniprot ID</b>         |
| <b>Immunogen</b>          |
| <b>Immunogen Region</b>   |
| <b>Specificity</b>        |
| <b>Immunogen Sequence</b> |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081