

**Anti-p120/Catenin antibody [MRQ-5] (STJ16100406)**  
STJ16100406**GENERAL INFORMATION**


|                          |                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                             |
| <b>Short Description</b> | Mouse monoclonal antibody anti-p120/Catenin is suitable for use in Immunohistochemistry research applications. |
| <b>Applications</b>      | IHC-P                                                                                                          |
| <b>Host/Source</b>       | Mouse                                                                                                          |
| <b>Reactivity</b>        | Human                                                                                                          |

**PRODUCT PROPERTIES**

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>Clone ID</b>            | MRQ-5                                                             |
| <b>Concentration</b>       |                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Purification</b>        |                                                                   |
| <b>Dilution Range</b>      | IHC 1:50 – 1:200                                                  |
| <b>Formulation</b>         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| <b>Isotype</b>             | IgG1                                                              |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt.                                      |

**TARGET INFORMATION**

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen Region**  
**Specificity**  
**Immunogen Sequence**



---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081