

Anti-Cytomegalovirus antibody [8B1.2, 1G5.2&2D4.2] (STJ16100451)

STJ16100451

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse mix antibody anti-Cytomegalovirus is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Mix
Clone ID	8B1.2, 1G5.2&2D4.2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:25 – 1:50
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081