

Anti-MadCAM-1 antibody [314G8] (STJ16100465)

STJ16100465

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MadCAM-1 is suitable for use in Immunohistochemistry, ELISA, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/ELISA/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	314G8
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:50 FC 1:50 IHC 1:50 ELISA 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	8174
Gene Symbol	MADCAM1
Uniprot ID	MADCA_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	