

Anti-Nectin-3 antibody [103-A1] (STJ16100466)

STJ16100466

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Nectin-3 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence and Immunoprecipitation research applications.
Applications	IHC-F/FC/ELISA/IF/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	103-A1
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	25945
Gene Symbol	NECTIN3
Uniprot ID	NECT3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	