

Anti-Actin antibody [647] (STJ16100467)

STJ16100467

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Actin is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/IHC-P/IF/WB
Host/Source	Mouse
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	647
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Antigen isolated from chicken skeletal muscle cells
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	