

Anti-Beta-Tubulin antibody [TU-06] (STJ16100474)

STJ16100474

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta-Tubulin is suitable for use in Western Blot, Immunohistochemistry and Immunocytochemistry research applications.
Applications	WB/IHC-P/ICC
Host/Source	Mouse
Reactivity	Arabidopsis/Mouse/Rat/Pig/Nicotiana/Paramecium/Chicken/Human/Fish/Tetrahymena

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TU-06
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Beta-subunits of porcine brain tubulin.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081