

**Anti-LAT antibody [LAT-01] (STJ16100483)**

STJ16100483

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-LAT is suitable for use in Western Blot, Immunoprecipitation and Flow Cytometry research applications.
<b>Applications</b>	WB/IP/FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	LAT-01
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">27040</a>
<b>Gene Symbol</b>	<a href="#">LAT</a>
<b>Uniprot ID</b>	<a href="#">LAT_HUMAN</a>
<b>Immunogen</b>	Bacterially produced recombinant polypeptide corresponding to the entire cytoplasmic domain of human LAT.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	